

Proline (Pro) Colorimetric Assay Kit

Catalog Number: AKAS0504

Price: [Online Inquiry](#)

Product Information	
Product Name	Proline (Pro) Colorimetric Assay Kit
Detection Method	Colorimetric method
Instruments	Microplate Reader
Sample Type	Plant Tissue And Honey Sample
Detection Principle	Proline can react with acidic-ninhydrin to form stable red compound. The maximum absorption peak of the compound is at 520 nm. And the concentration of Pro can be calculated by measuring the OD value at 520 nm
Other Reagents	Acetic acid
Storage	2-8°C