

## CLO(Clorprenaline) ELISA Kit

**Catalog Number:** FSK0053

**Price:** [Online Inquiry](#)

Product Information	
Product Name	CLO(Clorprenaline) ELISA Kit
Reactivity	Universal
Sample Type	Muscle, Serum, Urine
Analytes	Quantitative analysis of Clorprenaline.
Reaction Temperature	25°C
Assay Time	30 min, 15 min
LOD (Detection Limit)	Serum-0.25 ppb; Urine, Muscle-0.5 ppb
Sample Preparation	Preparation methods are included in the kit description.
Instruments	Spectrophotometer, Microplate, Auto-analyser
Detection Method	Absorbance
Wave length (nm)	450 nm/630 nm
Detection Range	See the manual for details
Storage	2-8°C, avoid freeze
Expiration Date	1 year