

Mitochondrial Extraction Kit

Catalog Number: AKAS0422

Price: [Online Inquiry](#)

Product Information	
Product Name	Mitochondrial Extraction Kit
Sample Type	Animal Tissue
Detection Principle	The principle of mitochondrial separation and purification: after rupturing cells mechanically, the debris and huge organelles are removed by low-speed differential centrifugation, then purified mitochondria will be obtained by high-speed differential centrifugation.
Other Reagents	Double distilled water, PBS(0.01 M , pH=7.4)
Storage	-20℃