

CD Human P-Selectin Assay

Catalog Number: AK0365

| Catalog Number 2 | Size | Price in USD |
|------------------|---------------|--------------------------------|
| AK0365-1 | 1 plate size | Online Inquiry |
| AK0365-2 | 5 plate size | Online Inquiry |
| AK0365-4 | 25 plate size | Online Inquiry |

Product Information

| | |
|--------------------|--|
| Product Name | CD Human P-Selectin Assay |
| Description | The CD platform offers validated components that can be combined to create high quality assays quickly and easily. CD Human P-Selectin Assay has the essential components for an assay that measures P-Selectin in a single-analyte format. For information on multiplexing P-Selectin with other analytes, please see CD Custom Immuno-Oncology Group 1 (human) Assays. Related Key Terms: CD62, GRMP, PSEL, CD62P, GMP140, LECAM3, PADGEM, SELP, inflammation, leukocyte recruitment, cancer |
| Applications | Oncology & Cancer |
| Analytes | P-Selectin |
| Species | Human |
| Instruments | SECTOR Imager 2400 , MESO QuickPlex SQ 120MM , MESO SECTOR S 600 , SECTOR Imager 6000 , MESO SECTOR S 600MM , MESO QuickPlex SQ 120 |
| Sensitivity | 0.88 pg/mL |
| Dynamic Range | 0.88 - 15,000 pg/mL |
| Sample Type | Normal Lysate, Plasma, Serum, Tumor Lysate |
| Gene ID | 6403 |
| UniProt ID | P16109 |
| Schedule B Code | 3002.15.0100 |
| Plate Type | 96-well |
| Usage | For Research Use Only. Not for use in diagnostic procedures. |
| Storage Conditions | Multi-Component |

Kit Contents

AK0365-1

| Product Description | Storage | Quantity per Kit |
|---------------------|---------|------------------|
|---------------------|---------|------------------|

| Kit Contents | | |
|--|---------|----------------------|
| Calibrator 29 | 2-8 °C | 1 vial(1 vial) |
| MSD GOLD 96 Small Spot SA Plate Pack, SECTOR Plate | 2-8 °C | 1 each(1 plate) |
| Diluent 100 | 2-8 °C | 1 bottle(50 mL) |
| Diluent 3 | ≤-10 °C | 1 bottle(8 mL) |
| Diluent 58, 10 mL | ≤-10 °C | 1 bottle(10 mL) |
| MSD GOLD Read Buffer B | RT | 1 bottle(18 mL) |
| CD Human P-Selectin Antibody Set (1 Plate Size) | 2-8 °C | 1 each(1 plate size) |

AK0365-2

| Product Description | Storage | Quantity per Kit |
|--|---------|----------------------|
| Diluent 100 | 2-8 °C | 1 bottle(50 mL) |
| MSD GOLD Read Buffer B | RT | 1 bottle(90 mL) |
| Diluent 58 | ≤-10 °C | 1 bottle(50 mL) |
| Diluent 3 | ≤-10 °C | 1 bottle(40 mL) |
| MSD GOLD 96 Small Spot SA Plate Pack, SECTOR Plate | 2-8 °C | 5 each(1 plate) |
| Calibrator 29 | 2-8 °C | 5 vial(1 vial) |
| CD Human P-Selectin Antibody Set (5 Plate Size) | 2-8 °C | 1 each(5 plate size) |

AK0365-4

| Product Description | Storage | Quantity per Kit |
|--|---------|----------------------|
| CD Human P-Selectin Antibody Set (5 Plate Size) | 2-8 °C | 5 each(5 plate size) |
| Calibrator 29 | 2-8 °C | 25 vial(1 vial) |
| MSD GOLD 96 Small Spot SA Plate Pack, SECTOR Plate | 2-8 °C | 25 each(1 plate) |
| Diluent 100 | 2-8 °C | 5 bottle(50 mL) |
| Diluent 3 | ≤-10 °C | 5 bottle(40 mL) |
| Diluent 58 | ≤-10 °C | 5 bottle(50 mL) |
| MSD GOLD Read Buffer B | RT | 5 bottle(90 mL) |