

Rat NT-proANP Kit

Catalog Number: AK0060

Catalog Number 2	Size	Price in USD
AK0060-1	1 plate	Online Inquiry
AK0060-2	5 plates	Online Inquiry
AK0060-4	25 plates	Online Inquiry

Product Information	
Product Name	Rat NT-proANP Kit
Applications	Cardiac Biomarkers , Cardiac Injury , Natriuretic Peptide , Toxicology , Cardiovascular
Analytes	NT-proANP
Species	Rat
Instruments	SECTOR Imager 2400 , MESO SECTOR S 600MM , SECTOR Imager 6000 , MESO QuickPlex SQ 120MM , MESO QuickPlex SQ 120 , MESO SECTOR S 600
Schedule B Code	3002.15.0100
Plate Type	96-well
Usage	For Research Use Only. Not for use in diagnostic procedures.
Storage Conditions	Multi-Component

Kit Contents

AK0060-1

Product Description	Storage	Quantity per Kit
Rat NT-proANP Calibrator	≤-70 °C	1 vial(20 µL)
SULFO-TAG Anti-rat NT-proANP Antibody	2-8 °C	1 vial(75 µL)
Rat NT-proANP Plate	2-8 °C	1 each(1 each)
Diluent 30	≤-10 °C	1 bottle(25 mL)
Read Buffer T4X	RT	1 bottle(50 mL)

AK0060-2

Product Description	Storage	Quantity per Kit
Rat NT-proANP Calibrator	≤-70 °C	5 vial(20 µL)
SULFO-TAG Anti-rat NT-proANP	2-8 °C	1 vial(375 µL)

Kit Contents

Antibody		
Rat NT-proANP Plate	2-8 °C	5 each(1 each)
Diluent 30	≤-10 °C	5 bottle(25 mL)
Read Buffer T4X	RT	1 bottle(50 mL)

AK0060-4

Product Description	Storage	Quantity per Kit
Rat NT-proANP Calibrator	≤-70 °C	25 vial(20 µL)
SULFO-TAG Anti-rat NT-proANP Antibody	2-8 °C	5 vial(375 µL)
Rat NT-proANP Plate	2-8 °C	25 each(1 each)
Diluent 30	≤-10 °C	25 bottle(25 mL)
Read Buffer T4X	RT	5 bottle(50 mL)