

## Glycerol GK Assay Kit

**Catalog Number:** CAK0050

**Price:** [Online Inquiry](#)

Product Information	
Product Name	Glycerol GK Assay Kit
Content	70 assays / 700 assays / 600 assays
Note for Content	If all volumes are halved, the number of manual tests per kit can be doubled.
Description	This kit is used for for the determination and analysis of glycerol in beverages, foodstuffs and other material.
Sample Type	Food products, cosmetics, pharmaceuticals and other materials, etc.
Analytes	Glycerol
Instruments	Spectrophotometer, Microplate, Auto-analyser
Detection Method	Absorbance
Wavelength (nm)	340
Signal Response	Increase
Linear Range	1.0 to 35 µg of glycerol per
Limit of Detection	0.37 mg/L
Assay Time	7 min
Method of Detection	Novel method
Advantages	Increased stability  Competitive price (cost per test)  Short reaction time Standard included
Shipping Temperature	Ambient
Expiration Date	> 2 years
Storage Temperature	See individual component labels
Signal Word	Not Applicable
Hazard Statements	Not Applicable
Precautionary Statements	Not Applicable