

Starch Damage Assay Kit

Catalog Number: CAK0008

Price: [Online Inquiry](#)

Product Information	
Product Name	Starch Damage Assay Kit
Content	200 assays per kit
Description	This kit is used for the determination of starch damage in wheat flour / cereal flours.
Sample Type	Cereal.
Analytes	Starch Damage
Instruments	Spectrophotometer
Detection Method	Absorbance
Wavelength (nm)	510
Signal Response	Increase
Linear Range	4 to 100 µg of D-glucose per
Limit of Detection	0.5 g/100 g
Assay Time	40 min
Method of Detection	AACC Method, ICC Standard and RACI Standard Method
Shipping Temperature	Ambient
Expiration Date	> 1 years
Storage Temperature	See individual component labels
Signal Word	Danger
Hazard Statements	H314, H315, H319, H334
Precautionary Statements	P260, P261, P264, P280, P284, P301+P330+P331, P302+P352, P303+P361+P353, P304+P340, P342+P311, P501