

Total Starch Assay Kit (AA/AMG)

Catalog Number: CAK0012

Price: [Online Inquiry](#)

Product Information	
Product Name	Total Starch Assay Kit (AA/AMG)
Content	50 assays / 100 assays
Description	This kit is used for the determination of total starch in cereal flours and food products.
Sample Type	Cereal, food products samples.
Analytes	Total Starch
Instruments	Spectrophotometer
Detection Method	Absorbance
Wavelength (nm)	510
Signal Response	Increase
Linear Range	4 to 100 µg of D-glucose per
Limit of Detection	0.18 g/100 g total starch “as is”
Assay Time	90 min
Method of Detection	AACC Method 76-13.01, AOAC Method 996.11, ICC Standard Method No. 168 and RACI Standard Method
Shipping Temperature	Ambient
Expiration Date	> 2 years
Storage Temperature	See individual component labels
Signal Word	Danger
Hazard Statements	H314, H315, H319, H334
Precautionary Statements	P260, P261, P264, P280, P284, P301+P330+P331, P302+P352, P303+P361+P353, P304+P340, P342+P311, P501