

Carrez Clarification Kit

Catalog Number: CAK0092

Price: Online Inquiry

| Product Information | |
|----------------------------|---|
| Product Name | Carrez Clarification Kit |
| Content | 40-1000 assays per kit |
| Note for Content | This method is not applicable to sample preparation prior to the determination of ascorbic acid (vitamin C), citric acid and urea/ammonia. |
| Description | This kit is used for the removal of interfering compounds such as proteins and fats from samples, allowing improved handling and more accurate determination of target analytes such as carbohydrates and alcohols. |
| Sample Type | Dairy products samples. |
| pH Range | 7.5-8 (recommended) |
| Assay Time | 15 min |
| Advantages | Allows easy sample clarification of complex matrices with high protein and fat content |
| Shipping Temperature | Ambient |
| Expiration Date | > 2 years |
| Storage Temperature | See individual component labels |
| Signal Word | Danger |
| Hazard Statements | H318, H410 |
| Precautionary Statements | P273, P280, P305+P351+P338, P310, P391 |