

Total Sulfite Assay Kit (Enzymatic)

Catalog Number: CAK0081

Price: [Online Inquiry](#)

Product Information	
Product Name	Total Sulfite Assay Kit (Enzymatic)
Content	50 assays / 500 assays / 588 assays
Description	This kit is used for the determination of total sulfite (as SO ₂) in wines, beverages, foodstuffs and other materials. A simple, highly specific and reliable, enzymatic method.
Sample Type	Wine, beer, soft drinks, jam, fruit and vegetables, pharmaceuticals, etc.
Analytes	Total Sulfite
Instruments	Spectrophotometer, Microplate, Auto-analyser
Detection Method	Absorbance
Wavelength (nm)	340
Signal Response	Decrease
Linear Range	1 to 50 µg of sulfite per
Limit of Detection	0.34 mg/L
Assay Time	30 min
Method of Detection	DIN, EN, MEBAK and NMKL Method
Advantages	Extended cofactors stability Very cost effective Simple format Standard included
Shipping Temperature	Ambient
Expiration Date	> 2 years
Storage Temperature	See individual component labels
Signal Word	Not Applicable
Hazard Statements	Not Applicable
Precautionary Statements	Not Applicable