

## Glucuronoyl Esterase Assay Kit

Price: [Online Inquiry](#)

| Product Information      |                                                                                                                                                          |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------|
| Catalog No.              | EAK0011                                                                                                                                                  |
| Product Name             | Glucuronoyl Esterase Assay Kit                                                                                                                           |
| Content                  | 100 assays per kit                                                                                                                                       |
| Sample Type              | Fermentation broths, industrial enzyme preparations and biofuels research.                                                                               |
| Analytes                 | Glucuronoyl Esterase                                                                                                                                     |
| Instruments              | Spectrophotometer, Auto-analyser                                                                                                                         |
| Detection Method         | Absorbance                                                                                                                                               |
| Wave length (nm)         | 400                                                                                                                                                      |
| Signal Response          | Increase                                                                                                                                                 |
| Assay Time               | 10 min                                                                                                                                                   |
| Reproducibility          | 0.03                                                                                                                                                     |
| Method of Detection      | Novel method                                                                                                                                             |
| Advantages               | Specific for glucuronoyl esterase<br>Suitable for use in crude enzyme preparations<br>Suitable for manual and auto-analyser formats<br>Standard included |
| Shipping Temperature     | Ambient                                                                                                                                                  |
| Expiration Date          | > 2 years                                                                                                                                                |
| Storage Temperature      | See individual component labels                                                                                                                          |
| Signal Word              | Danger                                                                                                                                                   |
| Hazard Statements        | H300, H334, H411                                                                                                                                         |
| Precautionary Statements | P261, P264, P270, P273, P284, P301+P310, P304+P340, P342+P311, P501                                                                                      |